




Professional
HOSPITAL FURNISHERS

Single Use Instruments



Professional Hospital Furnishers

19-A (N) Industrial Estate,
P.O.Box: 2095, Sialkot-51310, Pakistan.
Tel.: +92 52 325 6441, 355 9980
Fax: +92 52 325 6442


www.phfworld.com

sale@phfworld.com



Index

- 1-Probes & Suction Tubes
- 2-Scissors
- 3-Scalpel Handles & Knives
- 4-Tissue Forceps
- 5-Sponge Holding Forceps
- 6-Hemostatic Forceps
- 7-Needle Holders
- 8-Gynecological Instruments
- 9-Hooks & Retractors
- 10-Anesthesia Products
- 11-Eye Instruments
- 12-Dental Instruments
- 15-Supplementary Products



Jobson Horne Probe

Item code	Size in mm
SB-19.35	125
SB-19.38	180

Single Use Jobson Horne Probe, designed for the removal of earwax by a clinician.



Formby Cerumen Hook & Scoop 14cm

Item code	Size in mm
SAE-09.34	140

Single Use Formby Cerumen Hook & Scoops are used to remove foreign objects from the ear. Most commonly they are used for the removal of earwax.



Wagener Ear Hook, Fig.1

Item code	Size in mm
SAC-79.01	165

Single Use Hook used to remove foreign objects from the ear. The removal of earwax is the most common use of this instrument.



Zoellner Suction Tube 18cm

Item code	Size in mm
SAL-01.18	180

Single Use Zoellner Suction Tubes, used to remove fluids from the operating area during ENT procedures.

Silver Probe with Eye 13cm

Item code	Size in mm
SB-06.13	130

Single Use Silver Probe with Eye, designed to gauge the depth or direction of a sinus or cavity.





Iris Scissors

Item code	Size in mm	
SC-01.11	115	Straight
SC-02.11	115	Curved

Single Use Straight Iris Stitch Scissors, used to cut suture material or delicate tissue during minor surgical procedures.



Mayo Scissors

Item code	Size in mm	
SC-24.15	155	Straight
SC-24.17	170	Straight
SC-25.17	170	Curved
SC-25.23	230	Curved

Single Use Straight Mayo Scissors, used to cut and dissect tough tissue, or cut sutures.



Operating Scissors

Item code	Size in mm	
SC-14.13	130	Blunt/Blunt
SC-15.13	130	Sharp/Sharp
SC-16.13	130	Sharp/Blunt

Single Use Dressing Scissors, used as a general scissor to cut small and medium wound dressings.



Lister Bandage Scissors

Item code	Size in mm	
SL-12.14	145	
SL-12.18	180	

Single Use Lister Bandage Scissors, used to cut thick or heavy bandages and dressings.

Fergusson Abdominal Angled on Flat

Item code	Size in mm
SC-97.20	200



Single Use Currie Uterine Scissors designed to make an incision to the perineum and used in episiotomy procedures.

Strabismus Scissor

Item code	Size in mm	
SAB-196.42	115	Straight
SAB-196.43	115	Curved



Single Use Curved Strabismus Scissors with blunt/blunt blade ends, used specifically during the treatment of Strabismus (ocular misalignment) and other types of eye surgery.

Metzenbaum Scissors

Item code	Size in mm	
SC-36.14	140	Straight
SC-37.14	140	Curved



Single Use Metzenbaum Scissors, used for the dissecting and cutting of delicate tissue during minor surgical procedures.

Metzenbaum Scissors Fine

Item code	Size in mm
SC-42.14	140



Single Use Fine Straight Metzenbaum Scissors, used for the dissecting and cutting of delicate tissue during minor surgical procedures.



Metzenbaum Scissors

Item code	Size in mm	
SC-38.18	180	Straight
SC-39.18	180	Curved

Single Use Straight Metzenbaum Scissors, used for the dissecting and cutting of delicate tissue during minor surgical procedures.



Kilner Scissors

Item code	Size in mm	
SC-59.11	115	Straight
SC-59.13	115	Curved

Single Use Straight Kilner Scissors, used for the dissecting and cutting of delicate tissue during minor surgical procedures.



McIndoe Scissors

Item code	Size in mm	
SC-78.18	180	Straight
SC-79.18	180	Curved

Single Use Straight McIndoe Scissors, used for the dissecting and cutting of tissue during surgical procedures



Umbilical Scissors, Usa Pattern, 10.5cm

Item code	Size in mm
SP-32.11	105

Single Use Umbilical Cord Scissors, used to cut the umbilical cord during child birth.

Suture Scissors Hook Point 11.5cm

Item code	Size in mm
SS-20.12	115

Single Use Hook Point Suture Scissors, designed specifically for the cutting and removal of sutures from patients.



Tough Cut Scissors Non-Sterile 18.5cm

Item code	Size in mm
SL-08.18	180

Non-Sterile Tough Cut Scissors, used to cut heavy or thick bandages and for the removal of light casts.



Tenotomy Scissors

Item code	Size in mm
SL-08.18	114.3

Tenotomy Scissors, Curved Tip
Surgical Stainless Steel





Scalpel Handle BP No.3

Item code	Size in mm
SE-04.03	125

Single Use BP Scalpel Handle No. 3, designed to fit scalpel blade sizes 10, 10A, 11, 12 and 15.



Scalpel Handle BP No.9

Item code	Size
SE-10.84	#9

Single Use BP Scalpel Handle No. 9, designed to fit scalpel blade sizes 10-15.

Arthur Splinter Forceps 14cm

Item code	Size in mm
SF-88.24	140

Single Use Arthur Splinter Forceps, used for the removal of splinters and other foreign bodies, particularly from wounds to the hands or feet.



Martins Splinter Forceps 11.5cm

Item code	Size in mm
SF-91.11	115

Single Use Martins Splinter Forceps, used for the removal of splinters and other foreign bodies from a wound, or to hold and manoeuvre dressings or soft tissue during procedures.



Dental Forceps 23cm

Item code	Size in mm
SH-10.52	140

Single Use Blunt Vasectomy Forceps, used to grasp and hold the Vas during a vasectomy procedure.



McIndoe Forceps Non-Toothed 15cm

Item code	Size in mm
-----------	------------

SF-35.15	150
----------	-----

Single Use McIndoe Dissecting Forceps, used to grasp and hold tissue or foreign bodies in a surgical wound. Also used to hold wound dressings and packing.

McIndoe Forceps Toothed 15cm

Item code	Size in mm
-----------	------------

SF-35.85	150
----------	-----

Single Use McIndoe Toothed Dissecting Forceps – 1:2 Teeth, used to grasp and hold tissue in a surgical wound.

DeBakey Dissecting Forceps

Item code	Size in mm
-----------	------------

SF-43.15	160
SF-43.19	190

Single Use DeBakey Dissecting Forceps, used for surgical procedures to hold vascular vessels and tissue.

Hunter Splinter Forceps 11cm

Item code	Size in mm
-----------	------------

SF-93.11	110
----------	-----

Single Use Hunter Splinter Forceps, used for the removal of splinters and other foreign bodies from a wound, or to hold and manoeuvre dressings or soft tissue during procedures.

Adson Forceps Non-Toothed 12.5cm

Item code Size in mm

SF-15.12 125

Single Use Adson Tissue Forceps, used for fine surgical procedures to hold delicate or superficial tissues. Also used to tie sutures at the end of the procedure and hold dressings.



Adson Forceps Toothed 12.5cm

Item code Size in mm

SF-16.12 125

Single Use Adson Toothed Tissue Forceps – 1:2 Teeth, used for fine surgical procedures to hold delicate or superficial tissues.



Adson Micro Forceps Non-Toothed 12.5cm

Item code Size in mm

SF-17.12 125

Single Use Adson Micro Tissue Forceps, used for very fine surgical procedures to hold very delicate or superficial tissues or hold dressings.



Adson Micro Forceps Toothed 12.5cm

Item code Size in mm

SF-19.12 125

Single Use Adson Toothed Micro Tissue Forceps – 1:2 Teeth, used for very fine surgical procedures to hold very delicate or superficial tissues.





TOE Forceps 13cm

Item code	Size in mm
-----------	------------

SF-113.13	130
-----------	-----

Single Use Turn Over End (TOE) Dissecting Forceps, designed to grasp or hold soft tissues during procedures. Also used to remove debris from wounds and aid the application of dressings.



Treves Toothed Forceps, 1x2 Teeth 13cm Non-Toothed

Item code	Size in mm
-----------	------------

SF-23.13	130
----------	-----

Single Use Treves Toothed Dissecting Forceps – 1:2 Teeth, designed to grasp or hold soft tissues during procedures. Also used to aid the application of dressings.



Iris Forceps Non-Toothed / Curved 10cm

Item code	Size in mm
-----------	------------

SF-24.10	100
----------	-----

Single Use Iris Non-Toothed Dissecting Forceps, used for fine surgical procedures to hold delicate tissues during suturing. Also used to tie sutures at the end of the procedure and hold dressings.



Iris Forceps Toothed, 1x 2 teeth 10cm Bonney Forceps 1x 2

Item code	Size in mm
-----------	------------

SF-27.10	100
----------	-----

Single Use Iris Toothed Dissecting Forceps – 1:2 Teeth, used for fine surgical procedures to hold delicate tissues during suturing. Also used to tie sutures at the end of the procedure.

Ramsey Sponge Holder 18cm

Item code	Size in mm
SG-25.18	180

Single Use Ramsey Sponge Holder, suitable for holding swabs when 'prepping' a patient or used to grip and manipulate cervical tissue during gynaecological procedures



Ramsey Sponge Holder 25cm

Item code	Size in mm
SG-37.24	250

Single Use Ramsey Sponge Holder, suitable for holding swabs when 'prepping' a patient or used to grip and manipulate cervical tissue during gynaecological procedures



Lister Sinus Forceps 13cm

Item code	Size in mm
SG-01.13	130

A single use Lister Sinus Forceps used for the removing or inserting of dressings to the sinus.





Backhaus Towel Clip 10cm

Item code Size in mm

SI-25.10	100
SI-26.13	130

Single Use Backhaus Towel Clip, designed to secure towels and surgical draping during a medical procedure.



Spencer Wells Artery Forceps

Item code Size in mm

SH-54.13	130	Straight
SH-54.16	160	Straight
SH-54.18	180	Straight
SH-54.20	200	Straight
SH-54.23	230	Straight
SH-55.16	160	Curved
SH-55.18	180	Curved

Single Use Straight Spencer Wells Artery Forceps, used during surgery to compress the artery, clamp and seal small blood vessels or hold the artery out of the way.



Allis Tissue Forceps

Item code Size in mm Teeth

SN-02.15	150	4 x 5
SN-02.18	180	5 x 6
SN-02.19	190	5 x 6
SN-02.25	250	5 x 6

Single Use Allis Tissue Forceps, designed to grasp and hold organs, skin and fascia tissue, or secure operating materials such as tubes and drapes.



Pean Forceps

Item code Size in mm

SH-42.14	140	Straight
SH-43.14	140	Curved

Single Use Pean Forceps, used during surgery to compress the artery, clamp and seal small to medium size blood vessels or hold the artery out of the way.

Kocher Artery Forceps

Item code	Size in mm	
SH-66.14	140	Straight
SH-67.14	140	Curved

Single Use Kocher Artery Forceps, used during surgery to compress the artery, clamp and seal blood vessels or hold the artery out of the way.



Halsted Mosquito Artery Forceps

Item code	Size in mm	
SH-10.12	125	Straight
SH-11.12	125	Curved

Single Use Halsted Mosquito Artery Forceps, used as a ratcheted haemostat to secure and clamp bleeding vessels. Also used to grasp delicate tissues.



Crile Forceps

Item code	Size in mm	
SH-46.14	140	Straight
SH-47.14	140	Curved

Single Use Crile Forceps, used during surgery to compress the artery, clamp and seal small to medium size blood vessels or hold the artery out of the way.



Vasectomy Forceps Blunt 16cm

Item code	Size in mm	
SH-10.52	140	

Single Use Blunt Vasectomy Forceps, used to grasp and hold the Vas during a vasectomy procedure.





Olsen-Hegar Needle Holder 14cm

Item code	Size in mm
SJ-63.14	140

Single Use Olsen-Hegar Needle Holder, designed to hold fine suturing needles when closing wounds. Features cross hatched jaws to maintain grip on needles and cutting blades to cut through the suture thread after the procedure.



Derf Needle Holder 13cm

Item code	Size in mm
SJ-00.12	130

Single Use Derf Needle Holder, designed to hold fine suturing needles when closing wounds. Features cross hatched jaws to maintain grip on needles.



Halsey Needle Holder 13cm

Item code	Size in mm
SJ-14.13	130

Single Use Halsey Needle Holder, designed to hold very fine suturing needles when closing wounds. Features flat faced jaws to maintain grip on very fine needles.



Kilner Needle Holder 14cm

Item code	Size in mm
SJ-59.13	140

Single Use Kilner Needle Holder, designed to hold fine suturing needles when closing wounds. Features cross hatched jaws to maintain grip on needles.

Mayo Needle Holder 13cm

Item code	Size in mm
SJ-23.14	130
SJ-23.15	150

Single Use Mayo Needle Holder, designed to hold suturing needles when closing wounds. Features cross hatched jaws with positioning groove to maintain grip on the needle.





Cusco Vaginal Speculum

Item code	Size in mm	
SO-32.62	15 x 95	Virgin
SO-32.63	25 x 100	Small
SO-32.64	30 x 105	Medium
SO-32.65	30 x 120	Medium Long
SO-32.66	35 x 105	Large

Single use large size Cusco Vaginal Speculum, used to retract the vaginal wall during gynaecological examinations or procedures, such as vaginal surgery.



Wintertons Vaginal Speculum

Item code	Size in mm
SO-34.95	27 x 130

Single Use Wintertons Vaginal Speculum used to retract the vaginal wall during gynaecological examinations or procedures, such as vaginal surgery.



Simms Vaginal Speculum

Item code	Size in mm	
SO-21.10	70 x 25mm, 75 x 30mm	Fig. 1
SO-21.20	75 x 30mm, 80 x 35mm	Fig. 2
SO-21.30	80 x 35mm, 90 x 40mm	Fig. 3

Single Use Sims Vaginal Speculum used to retract the vaginal wall during gynaecological examinations or procedures. This is a single blade speculum designed for use when patients are examined in the Sim's (or lateral) position.

Hegar Dilator

Item code	Size in mm
SO-40.01	1 + 2
SO-40.03	3 + 4
SO-40.05	5 + 6


Single Use Hegar Dilator used to dilate the cervix.



IUD Hook 31cm

Item code	Size in mm
SP-06.61	Ø26 x 2


Single Use IUD Hook designed to probe for and grasp the intrauterine device during its removal.



Currie Uterine Scissors 20cm

Item code	Size in mm
SC-97.20	200

Single Use Currie Uterine Scissors designed to make an incision to the perineum and used in episiotomy procedures.



Simms Uterine Sound 32cm

Graduation in cm

Item code	Size in mm
SO-49.31	320

Single Use Currie Uterine Scissors designed to make an incision to the perineum and used in episiotomy procedures.





KOVOURKIAN Baby Tischler Cervical Punch Forceps 23cm

Item code	Size in mm
-----------	------------

SO-134.25	230 shaft
-----------	-----------

Single Use Baby Tischler Cervical Punch Biopsy Forceps, designed to help diagnose any irregularities of the cervix and the sample the endocervical canal with minimal patient discomfort and bleeding.



Schumacher Cervical Punch Forceps 23cm

Item code	Size in mm
-----------	------------

SO-128.01	230
-----------	-----

Single Use Schumacher Cervical Punch Biopsy Forceps used to take a section of tissue from the cervix as a sample during gynaecological procedures.



Bonney Uterine Polyp Forceps 24cm

Item code	Size in mm
-----------	------------

SG-26.34	240
----------	-----

Single Use Bonney Uterine Polyp Forceps used for the retrieval of polyps from the endometrial cavity of the uterus, or stones from the gallbladder or duct.



Umbilical Scissors, Usa Pattern, 10.5cm

Item code	Size in mm
-----------	------------

SP-32.11	105
----------	-----

Single Use Umbilical Cord Scissors, used to cut the umbilical cord during child birth.

Vulsellum Forceps 23cm 1x1Teeth

Item code	Size in mm
SO-63.26	230


Single Use Luer Vulsellum Forceps, suitable for grasping the uterus during gynaecological procedures. Alternatively referred to as Tenaculum Forceps.



Teale Vulsellum Forceps 23cm

Item code	Size in mm
SO-98.33	230

Single Use Teale Vulsellum Forceps, suitable for grasping the uterus during gynaecological procedures. Alternatively referred to as Tenaculum Forceps.





Gillies Skin Hook 18cm

Item code	Size in mm
-----------	------------

SY-12.20	180
----------	-----

Single Use Gillies Skin Hook, used to retract the skin during surgical procedures such as the dissection of lesions.



Magills Forceps

Item code	Size in mm	
SAD-05.19	190	Child
SAD-05.24	240	Adult

Magills Forceps, designed to hold cotton balls to swab the vocal cords before intubation. The 25cm version is primarily designed for use with adults. This is a single use instrument,





Dental Syringe 2.2mm Cartridge

Item code	Size
20.0703S	2.2CC cartridge

Single Use Dental Syringe for 2.2mm cartridges. Designed for administering an anaesthetic during dental or a variety of other procedures.





Hartman Crocodile Forceps 22cm

Item code Size in mm

SAC-103.47 220 Shaft

Single Use Hartmann Crocodile Forceps, used as grasping forceps to hold dressings and swabs, or to remove foreign bodies.



Hartman Crocodile Forceps 14.5cm Extra Fine

Item code Size in mm

SAC-103.48 145

Single Use Extra Fine Hartmann Crocodile Forceps, used for the insertion or removal of grommets and t-tubes, or to remove foreign bodies.



Trousseau Bowlby Tracheal Dilating Forceps 14cm

Item code Size in mm

SAP-36.14 140

Single Use Trousseau Bowlby Tracheal Dilating Forceps used within ENT procedures to dilate tracheostomy tracks.



Tilley Aural Dressing Forceps 14cm

Item code Size in mm

SAC-105.64 140

Tilley Aural Dressing Forceps are designed for the removal or insertion of nasal dressings or packs.



Thudicum Nasal Speculum

Item code	Size
SAM-17.00	#0
SAM-17.01	#1
SAM-17.02	#2
SAM-17.03	#3
SAM-17.04	#4
SAM-17.05	#5
SAM-17.06	#6
SAM-17.07	#7

Thudicum Nasal Speculum designed to open and expand the nasal cavity.



Volkman Spoon Double Ended

Item code	Size in mm
SW-24.13	130
SW-24.14	140
SW-24.17	170

Single Use Double Ended Volkman Spoon, designed to scrape away carious bone or other diseased tissue from wounds.



Volkman Spoon Double Ended 22cm

Item code	Size in mm
SW-25.01	220

Single Use Double Ended Volkman Spoon, designed to scrape away carious bone or other diseased tissue from wounds.

